


Introduction to XEDS in the AEM



Nestor J. Zaluzec

zaluzec@microscopy.com

zaluzec@aaem.amc.anl.gov

Brief Review of X-ray Generation

Instrumentation: Detector Systems


Instrumentation: EM Systems



Data Analysis and Quantification:

Additional Topics

Brief Review of X-ray Generation


Electron Excitation of Inner Shell & Continuum Processes
Characteristic and Bremsstrahlung Emission
Spectral Shapes
Notation of Lines






Schematic Diagram Illustrating Sources of Inelastic Scattering Signals

Experimental XEDS, XPS, and EELS data from the Copper L shell.
Note the differences in energy resolution, and spectral features.



Nomenclature for Principle X-ray Emission Lines



$$\text{Characteristic X-ray Line Energy} = E_{\text{final}} - E_{\text{initial}}$$

Recall that for each atom every shell has a unique energy level determined by the atomic configuration for that element.

X-ray line energies are unique.

Nomenclature for X-ray Lines


X-ray Transition Selection Rules: (Principle Quantum Numbers)

Shell	n	l	j	Rule
K	1	0	1/2	
L	2	0,1	1/2, 3/2	$\Delta n > 0$
M	3	0,1,2	1/2, 3/2, 5/2	$\Delta l = +1, -1$
N	4	0,1,2,3	1/2, 3/2, 5/2, 7/2	$\Delta j = +1, 0, -1$

Relative Intensities of Major X-ray Lines

$K_{\alpha 1} = 100$	$L_{\alpha 1} = 100$	$M_{\alpha 1,2} = 100$
$K_{\alpha 2} = 50$	$L_{\alpha 2} = 50$	$M_{\beta} = 60$
$K_{\beta 1} = 15-30$	$L_{\beta 1} = 50$	
$K_{\beta 2} = 1-10$	$L_{\beta 2} = 20$	
$K_{\beta 3} = 6-15$	$L_{\beta 3} = 1-6$	
	$L_{\beta 4} = 3-5$	
	$L_{\gamma 1} = 1-10$	
	$L_{\gamma 3} = 0.5-2$	
	$L_{\eta} = 1$	
	$L_{\zeta} = 1-3$	

Characteristic X-Ray Spectrum Illustrating KLM lines





TEM Specimen: Y₁Ba₂Cu₃O_{6.9} Superconductor - 120 kV - UTW Detector

Note:

As Z increases the Kth shell line energy increases.
If K-shell is excited then all shells are excited (Y, Cu, Ba)
but may not be detected.
Severe spectral overlap may occur for low energy lines.

Electron Excitation of Continuum Processes






Energy Range - Continuous Distribution

Maximum = Incident Electron Energy (Least Frequent)

Minimum = $E_{\text{plasmon}} \sim 15\text{-}30\text{ eV}$ (Most Frequent)

Spectral Distribution will reflect this range, modified by detector response function

Electron Excitation of Continuum (Background) Intensity



Notes:

Spectral background will be influenced by:

- 1.) Specimen composition
- 2.) Detector efficiency
- 3.) TEM generated artifacts

Instrumentation: Detector Systems

Wavelength Dispersive Spectrometers (WDS)
Energy Dispersive Spectrometers (EDS)
 Si(Li) Detectors
 HPGe Detectors
Spectral Artifacts of the EDS System
Detector Efficiency Functions
Light Element Detectors
Multichannel Analyzers

Energy Dispersive Spectrometers: (Solid State Detector)

Operates on Energy Deposition Principle

Simple, Nearly Operator Independent
Large Solid Angles (0.05-0.3 sr)
Virtually Specimen Position Independent
No Moving Parts
Parallel Detection
Quantification by Standardless or Standards Methods

Poor Energy Resolution (~ 130 eV)
 ** Superconducting Systems (~ 20 eV)
Poor Peak/Background Ratios (100:1)
Detection Efficiency Depends upon X-ray Energy




FIGURE 10. Operating schematic of a Si(Li) detector system.

Wavelength Dispersive Spectrometers : (Diffractometer)

Operates using Diffraction Principles (Bragg's Law)

Excellent Energy Resolution (~ 5 eV)
High Peak/Background Ratios (10000:1)
Good Detection Efficiency for All X-rays
High Counting Rates
Good Light Element Capabilities

Complex Mechanical Devices, Operator Intensive
Specimen Height dependant focus
Moving Components in the AEM
Limited Solid Angles (<0.01 sr)
Serial Detection
Quantification Requires Standards






FIGURE 1. Basic components of a crystal diffraction spectrometer system.

Comparison of EDS and WDS Spectrometers


<u>Parameter</u>	<u>Wavelength Dispersive</u>	<u>Energy Dispersive</u>
Construction	Mechanical Device moving components	Solid State no moving parts
Energy Resolution	5 eV	130 eV
Efficiency	≤ 30 %	100 % (3-15keV)
Input Count Rate	30-50 K cps	10 K cps
Peak/Background*	10000	100
Atomic Number Range	$Z \geq 4$ (Be)	$Z \geq 11$ (Na) $Z \geq 5$ (B)
Number of Elements	1 per Detector	All in Energy Range
Solid Angle	0.001-0.01 sr	0.02-0.3 sr
Collection Time	Tens of Minutes	Minutes
Beam Current	High Stability Required	Low Stability Required
Detector Stability	Good Short Term	Excellent
Spectral Artifacts	Negligible	Important
Operation	Skilled (?)	Novice

* Values depend on definition, specimen, and operating conditions

Comparison of EDS and WDS Spectra




Operation of a Crystal Spectrometer
Using Braggs' Law $n\lambda = 2d \sin(\theta)$



Installation of a Crystal Spectrometer in a TEM - EMMA-4 System





• P-type Silicon high conductivity due to impurities (usually Boron); Lithium acts as a compensating dopant neutralizing the Si giving it a high resistivity.


• Radiation deposits energy in the Si(Li) lattice & creates free electron-hole pairs in the crystal @ 1 electron-hole pair/3.8 eV of deposited energy @ 77K.

• Intrinsic semiconducting Si allows both electrons & holes to become mobile under application of a potential bias across the crystal

Properties of Intrinsic Silicon

Attaching HV electrodes to the two surfaces the Si(Li) crystal will act similar to a capacitor with free charges developing on the electrical contacts. Charge developed in the crystal is $N = E / \dots$ (E= x-ray Energy, = 3.8 eV/e-h pair) i.e. \Rightarrow 10 kV X-ray produces ~ 2630 electrons = 4.2×10^{-16} Coulombs.

Solid State Detector Construction




Relative Detection Efficiency

Solid State Detectors: Si(Li) or Intrinsic (High Purity) Ge
Using a simple absorption model define the relative detector efficiency (E) by the following procedure:

$$I_o \quad \text{---} \quad x \quad \text{---} \quad I_T = I_o \exp(-\mu x) = I_o \exp\left(-\left[\frac{\mu}{\rho}\right] \rho x\right)$$

Components of Si(Li) Detector



Relative Detection Efficiency

Solid State Detectors: Si(Li) or Intrinsic (High Purity) Ge
Using a simple absorption model define the relative detector efficiency (ϵ) by the following procedure:


$$I_T = I_o \exp(-\mu x) = I_o \exp\left(-\left[\frac{\mu}{\rho}\right] \rho x\right)$$

$$\epsilon(E) = \frac{I_A * I_T}{I_o} = \exp\left(\sum_i^{HC/Be/Au/DL} -\left(\frac{\mu(E)}{\rho}\right)_i * \rho_i * t_i\right) * \{1 - \exp\left(-\left(\frac{\mu(E)}{\rho}\right)_j * \rho_j * t_j\right)\}$$

<-- Absorption --> <-- Transmission -->

$\mu(E)$ = mass absorption coefficient for Energy E; ρ = density; t = layer thickness


Calculated Si(Li) Detector Efficiency by Active Layer Thickness & Window Type



Relative Transmission Efficiency


$$\varepsilon(E) = \frac{I_r}{I_o} = \exp\left(-\sum_i^{\text{Window}} \left(\frac{\mu(E)}{\rho}\right)_i \rho_i t_i\right)$$

$\mu(E)$ = mass absorption coefficient for Energy E; ρ = density; t = layer thickness



Note the Variation in transmission characteristics by Window Type. Not all UltraThin Windows are Equivalent!!! For example Detection of Nitrogen using a Diamond window is virtually impossible.

Windowless vs. Conventional Detectors **Comparison of XEDS measurement on NiO** **using a Windowless versus Beryllium Window detector**



Note the enhanced detection efficiency below 1 keV for the WL detector. Both spectra are normalized to unity at the Ni K α Line (7.48 keV)

Windowless vs. Conventional Detectors


K Shell Spectra using Windowless Detector

Boron -> Silicon

L Shell Spectra Using Windowless Detector


Titanium -> Zinc

Note Potential Overlaps with K shell Lines




Comparison Light Element Spectroscopy Resolution XEDS vs EELS



Comparison of WL XEDS Detector and EELS spectra
taken from the same NiO specimen




Note the enhanced spectral information in the EELS data. Vertical scale is arbitrary and chosen for clarity of presentation.

Comparison Light Element Spectroscopy Resolution XEDS vs EELS






Resolution will also vary with
Microphonic & Electronic Noise, and Counting
Rate!




WL & UTW detectors are particularly sensitive to low energy noise and microphonics. Observe the changes in the spectra

Resolution Loss with Count Rate




Incomplete Charge Collection in Solid State Detectors




$$N_{\text{peak}} = N_0 \exp(-\mu d)$$

$$N_{\text{tail}} = N_0 (1 - \exp(-\mu d))$$

$$\frac{N_{\text{tail}}}{N_{\text{peak}}} = \frac{(1 - \exp(-\mu d))}{\exp(-\mu d)}$$




Operation of the Multi-Channel Analyzer



Adapted from Gedcke, 1972


MCA/ADC Considerations Detector Dead Time




Instrumentation: AEM Systems

The AEM as a system
Spectral Artifacts in the AEM
 Uncollimated Radiation
 Systems Peaks
 Artifacts at High Electron Energy
Specimen Contamination & Preparation
Optimizing Experimental Conditions


Spectral Artifacts in the AEM Uncollimated Radiation: The Hole Count




Spectral Artifacts in the AEM Uncollimated Radiation Solutions



Uncollimated Radiation: The Hole Count Effects of Thickness & Composition of Variable C2 Aperture




Spectral Artifacts in the AEM




Specimen: 200 Å Molybdenum Film on Holey Carbon supported on a Stnd Mo Aperture with 200 μm hole. Expt. Conditions: 120 kV, Specimen tilted toward Si(Li) 35 degrees.

Spectral Artifacts in the AEM


Hole Count Effects: Modified C_1 and C_2 Apertures



Hole Count Effects: Modified C_1 and C_2 & Non-Beam Defining Apertures



AEM/XEDS Hole Count Test Specimen





1989 Higgs Meeting XEDS/AEM Performance Panel Members


Bentley - Oak Ridge Nat. Lab
Cazaux - UFR Sciences
Craven - Univ. of Glasgow
Glas-Lab. de Bagneux
Hren-NC State Univ.


Kruit-Delft Univ.
Malisl-PMRL/CANMET
Rez-ASU
Williams-Lehigh Univ.
Zaluzec-Argonne Nat. Lab.

Spectral Artifacts in the AEM





Examples of X-ray Spectra from the AEM / XEDS Cr/Mo test specimen
a) beam on the specimen, b.) beam in the "hole".






**Spectral Artifacts in the AEM:
Electrons Entering the
Detector
50 kV**

**Note the Effects of 50 keV
Electrons Entering the Detector
on Background**



**Spectral Artifacts in the AEM:
Electrons Entering the
Detector
300 kV**

**Note the Effects of 300 keV
Electrons Entering the
Detector on Background**



Optimizing Experimental Conditions

Choice of X-ray Line

K- series
L- series
M- series

Detector/Specimen Geometry

Elevation Angle
Solid Angle

Detector Collimation

Choice of Accelerating Voltage

Relative Intensity
Peak/ Background
Systems Peaks/Uncollimated Radiation


Choice of Electron Source

Spatial Resolution
Tungsten Hairpin
LaB₆
Field Emission

Radiative Partition Function (Γ) Governs the Relative Intensities
Nominal Values (Varies slowly with Atomic Number)


K Shell	L Shell	M Shell
$K_{\alpha 1} = 100$	$L_{\alpha 1} = 100$	$M_{\alpha 1,2} = 100$
$K_{\alpha 2} = 50$	$L_{\alpha 2} = 50$	$M_{\beta} = 60$
$K_{\beta 1} = 15-30$	$L_{\beta 1} = 50$	
$K_{\beta 2} = 1-10$	$L_{\beta 2} = 20$	
$K_{\beta 3} = 6-15$	$L_{\beta 3} = 1-6$	
	$L_{\beta 4} = 3-5$	
	$L_{\gamma 1} = 1-10$	
	$L_{\gamma 3} = 0.5-2$	
	$L_{\eta} = 1$	
	$L_{\chi} = 1-3$	

Detector/Specimen Geometry



Designation	Elevation Angle θ_E	Azimuthal Angle θ_A	Manufacturer
Low	0°	45°	JEOL
	0°	90°	JEOL, Philips, VG
Intermediate	15-30°	90°	Philips, JEOL, Hitachi, VG
High	68-72°	0°	Hitachi, JEOL

Detector/Specimen Geometry




Detector/Specimen Geometry

Collection Solid Angle


$$\omega = A/R^2 = 0.3 - 0.001 \text{ sr}$$

A = Detector Active Area
 = 10-30 mm²
 R = Crystal to Specimen Distance
 = 10 - 50 mm


Subtending Solid Angle



Detection of System Peaks Effects of the Collimator & Stage



Detection & Removal of System Peaks




Removal of Stage System Peaks by use of Beryllium Gimbals
 Ge specimen 10,000 in Ge K' peak in both spectra
 Left Standard Single Tilt Cu Stage, Right Be Gimbal DT Stage

Reactive Gas Plasma Processing Applications to Analytical Electron Microscopy

Example:

- The figure at the right shows the results of contamination formed when a 300 kV probe is focussed on the surface of a freshly electropolished 304 SS TEM specimen.
- The dark deposits mainly consist of hydrocarbons which diffuse across the surface of the specimen to the immediate vicinity of the electron probe. The amount of the contamination is a function of the time spent at each location. Here the time was varied from 15 - 300 seconds.




Comparison Results on Electropolished 304 SS


• Untreated Specimen

• After 5 minutes Argon Processing


• After 5 minutes of additional Oxygen Processing



Specimen Contamination Effects



Electropolishing leaves residual Cl on surface



Ion Milling leaves redeposited Fe, Ni, Cu from SS Holder & Cu Washer

X-ray Production = Cross-section * electrons

$$\eta = \text{Probe current} = \frac{(\pi d_0 \alpha_0)^2 \mathcal{B}}{4}$$

$$\mathcal{B} = \text{Brightness} = \left(\frac{I_c}{\pi k T} \right) e V_r$$

$$V_r = \text{Relativistic Voltage} = V_0 \left(1 + \frac{e V_0}{2 m_0 c^2} \right)$$

$$V_0 (1 + 9.785 \times 10^{-7} V_0)$$

Non-Relativistic Cross-section Model



$$Q_K = \frac{a_K * b_K * \left\{ \ln \left(c_K * \frac{E_0}{E_c} \right) \right\}}{E_0 E_c}$$

Relativistic Cross-section Model

$$Q_K = \frac{a_K * b_K * \left\{ \ln \left(c_K * \frac{T_0}{E_c} \right) - \ln(1 - \beta^2) - \beta^2 \right\}}{T_0 * E_c}$$


$$T_0 = \frac{1}{2} m_0 c^2 \beta^2, \text{ where } \beta = v/c$$

Maximizing
The Signal



Effects on Intensity with Accelerating Voltage for constant Probe size parameters


Relative Number of K-shell Ionizations




$$I = Q * \eta$$



For identical probe diameters one has higher x-ray production at higher voltage due to the increase in beam current.
Alternatively, one can achieve the same statistical intensity for smaller probes at higher

Experimental Variation of Beam Current and X-ray Intensity with Voltage



Ni Ka X-ray Intensity





Al K intensity/nA/Å as a function of accelerating voltage for different specimen thicknesses.

Minimum Detectable Mass

$$\text{MDM} \sim \frac{k}{P_x I_o T} = \frac{k^*}{P_x J_o d_o^2 T} -$$

Minimum Mass Fraction

$$\text{MMF} \sim \frac{k}{\sqrt{[P_x (\frac{P}{B})_x] I_o T}} = \frac{k^*}{\sqrt{[P_x (\frac{P}{B})_x] J_o d_o^2 T}}$$

$k, k^* = \text{Constants}$

$P_x = \text{Characteristic Signal}$
from element X


$(P/B)_x = \text{Peak to Background ratio for element X}$

$I_o = \text{Incident electron flux}$

$J_o = \text{Incident electron current density}$



$d_o = \text{Probe diameter}$

$T = \text{Analysis time}$



	<u>Current</u>	<u>Future</u>
<u>XEDS</u>		
Sensitivity	MMF~ 10^{-3} MDM~ 10^{-20}	~ 10^{-4} ~ 10^{-22}
Quantification	2-10% same	
Instrumentation	Si(Li),HPGe	Si(Li),HPGe, UTW, WL
Limitations	Radiation Damage	
<u>EELS</u>		
Sensitivity	MMF~ 10^{-4} MDM~ 10^{-22}	~ 10^{-5} ~ 10^{-24}
Quantification	10-20% same	
Instrumentation	Serial/Parallel	
Limitations	Radiation Damage	
Future Directions	Electronic/Atomic Structure	

Spatial Resolution / Beam Spreading Monte Carlo Calculations



Monte Carlo Calculations of B (Newbury & Myklebust - 1979)

Element	Z	Thickness			
		10nm	50nm	100nm	500nm
Carbon	6	0.22	1.9	4.1	33.0
Aluminium	13	0.41	3.0	7.6	66.4
Copper	29	0.78	5.8	17.5	244.0
Gold	79	1.71	15.0	52.2	1725.0

Analytic Formulation (Elastic Scattering - Goldstein et al 1977)

$$B = 625 \frac{Z}{E_0} \sqrt{\frac{\rho}{A}} t^{3/2}$$

B = Beam Broadening [cm]


Z = Atomic Number

E_0 = Accelerating Voltage [kV]

ρ = Density [gms/cm³]

A = Atomic Weight

t = Thickness [cm]




Element	Z	10nm	50nm	100nm	500nm
Carbon	6	0.16	1.8	5.13	57.4
Aluminium	13	0.26	1.9	8.12	90.9
Copper	29	0.68	7.6	21.4	*
Gold	79	15.5	17.3	*	*

*model invalid at higher kV and/or high scattering angles

Comparison of Electron Sources

Type		Brightness β/V_0 A/cm ² /sr/eV	Source Size (μ m)	Energy Spread (eV)	Noise	Stability	Cohereny	Vacuum (Torr)
Thermionic	Hairpin	1	50	2-3	Low	Good	Low	<10 ⁻⁴
	Pointed	5	10			Fair	Moderate	<10 ⁻⁵
LaB ₆	Poly Crystal	10-30	10	~1	Low	Good	Moderate	<10 ⁻⁶
	Single Crystal	20-50	5					
Field Emission	Thermal Assist	100-500	<100Å	~.5	Fair	Moderate	High	<10 ⁻⁸
	Cold	100-1000		<.25		Fair		<10 ⁻¹⁰




Data Analysis and Quantification:

Spectral Processing
Thin Film Quantification Methods
Specimen Thickness Effects:
Absorption
Fluorescence

Spectral Processing : XEDS

Spectrum = Characteristic Peaks + Background



Data Reduction

Simple: Linear Background Fit & Integration

Curve Fitting: Non-Linear Background & Profile Matching

Frequency (Digital) Filtering: Background Suppression & Reference Spectra Fitting


Deconvolution: Fourier Method for Resolution Enhancement

Background Modeling

Simple - Linear and/or Polynomial Interpolation

Modeling - Parametric Fits of Analytic Expressions
Phenomenological Expressions
Modified Bethe Heitler Model
Digital Filtering - Mathematical Suppression

Spectral Processing : XEDS Simple Data Reduction



Note: Must use peak integrals (I) and not peak amplitudes (A)

Recall that for a Gaussian Peak

$$I = \int_{-\infty}^{+\infty} A \exp\left(\frac{-x^2}{2\sigma^2}\right) dx = \sqrt{2\pi} \sigma A$$


Hence for the ratio of Intensities

$$\frac{I_1}{I_2} = \frac{\sigma_1 A_1}{\sigma_2 A_2} \approx \frac{A_1}{A_2}$$


Spectral Processing : XEDS

Background Modeling : Power Law/Parametric Fits

$$\text{Bgnd} = \epsilon^* \left(A \left(\frac{E-E_0}{E} \right)^2 + B \left(\frac{E-E_0}{E} \right) + C \right)$$



Spectral Processing : XEDS Curve Fitting : Linear Modeling



Let

$$G_{ij} = \exp \left(\frac{-(E_j - E_i)^2}{2\sigma_{Ei}^2} \right)$$

then


$$\begin{aligned} Y_1 &= A_1 * G_{11} + A_2 * G_{21} + A_3 * G_{31} \\ Y_2 &= A_1 * G_{12} + A_2 * G_{22} + A_3 * G_{32} \\ Y_3 &= A_1 * G_{13} + A_2 * G_{23} + A_3 * G_{33} \end{aligned}$$

Using simple matrix algebra solve for A

- Fast and simple procedure
- Presumes operator knows all elements present

* System must be calibrated
(i.e. E_c and E_e must be accurately known)

Spectral Processing : XEDS Curve Fitting : Linear Modeling



Spectral Processing : XEDS
Curve Fitting : Non- Linear Modeling

Calculate and Minimize χ^2

$$\chi^2 = \sum_{k=1}^M \left(\sum_{i=1}^N \frac{(y_i - Y_k)^2}{y_i^2} \right) \quad \text{with } Y_k = A_k \exp \left(\frac{-(E_i - E_k)^2}{2\sigma^2(E_k)} \right)$$

and y_i = raw data

χ^2 is minimized by searching (E, σ , A) parameter space

Pattern Search

- Sequential method: mechanical iteration of each parameter until a local minimum is found.
- Simplex method: simultaneous variation of all parameters.

Gradient Search

- Evaluate

Multiple Least Squares with Derivative Reference

- Linearization of the above problem.
 $E \rightarrow E + \delta E \quad \sigma \rightarrow \sigma + \delta \sigma$

$$I'(E) = I(E) - \delta E \frac{\delta I(E)}{\delta E} + \left(\sigma \delta \sigma + \frac{(\delta E)^2}{2} \right) \frac{\delta^2 I(E)}{\delta E^2}$$





Figure 2. (a,b) plane and response surface showing a single, arbitrary, simplex.

Simplex search of Parameter Space after Fiori, NBS Special Pub. 604, 1979


Multiple Least Squares with Derivative Reference Minimizes Errors due to Calibration Drift




McMillian et al -
MAS 1984

Spectral Processing : XEDS Digital Filtering


Background Suppression by Mathematical modeling
- Replace Data by new spectra formed by the
following linear operation.




Operator independent
Introduces severe spectral distortion



Spectral Processing : XEDS Digital Filtering



X-ray Production




Quantitative Analysis using XEDS

For a thin specimen

$$I_A^{K\alpha} = \left\{ \sigma_A(E,Z) \Gamma_A \omega_A \right\} C_A \left\{ \frac{N_o \rho}{W_A} \right\} \left\{ \eta_o t \right\} \left\{ \epsilon_A \Omega \right\}$$

I_A	=	Measured x-ray intensity per unit area
	=	K^{th} -shell ionization cross-section
	=	K^{th} -shell fluorescence yield
	=	K^{th} -shell radiative partition function
W	=	Atomic Weight
N_o	=	Avagadro's number
	=	Density
C	=	Composition (At %)
ϕ	=	Incident electron flux
t	=	Specimen thickness
	=	Detector efficiency
	=	Detector solid angle

Ionization Cross-Section



For K Shells


$$Q_K = \frac{a_K * b_K * \left\{ \ln \left(c_K * \frac{T_o}{E_c} \right) - \ln(1-\beta^2) - \beta^2 \right\}}{T_o * E_c}$$

$$T_o = 1/2 m_o c^2 \beta^2, E_c = \text{Shell Excitation Energy}, \beta = v/c$$

For L Shells

$$Q_L = \frac{a_L * b_L * \left\{ \ln \left(c_L * \frac{T_o}{E_c} \right) - \ln(1-\beta^2) - \beta^2 \right\}}{T_o * E_c}$$

X-ray Fluorescence Yield has Systematic Variation With Atomic Number



K shell

K vs L shell

Relative Intensities of Major X-ray Lines

$$I_K = I_{K\alpha} + I_{K\beta}$$

$$I_L = I_{L\alpha} + I_{L\beta} + I_{L\gamma} + I_{L\eta} + I_{L\zeta}$$

$$I_M = I_{M\alpha} + I_{M\beta}$$

Line	Initial	Final	Line	Initial	Final	Line	Initial	Final
$K_{\alpha 1} \Rightarrow$	L ₃	K	$L_{\alpha 1} \Rightarrow$	M ₅	L ₃	$M_{\alpha 1} \Rightarrow$	N ₇	M ₅
$K_{\alpha 2} \Rightarrow$	L ₂	K	$L_{\alpha 2} \Rightarrow$	M ₄	L ₃	$M_{\alpha 2} \Rightarrow$	N ₆	M ₅
$K_{\beta 1} \Rightarrow$	M ₃	K	$L_{\beta 1} \Rightarrow$	M ₄	L ₂	$M_{\beta} \Rightarrow$	N ₆	M ₄
$K_{\beta 2} \Rightarrow$	N ₂	K	$L_{\beta 2} \Rightarrow$	N ₅	L ₃	$M_{\gamma} \Rightarrow$	N ₅	M ₃
$K_{\beta 3} \Rightarrow$	M ₂	K	$L_{\beta 3} \Rightarrow$	M ₃	L ₁			
			$L_{\beta 4} \Rightarrow$	M ₂	L ₁			
			$L_{\gamma 1} \Rightarrow$	N ₄	L ₂			
			$L_{\gamma 3} \Rightarrow$	N ₃	L ₁			
			$L_{\eta} \Rightarrow$	M ₁	L ₂			
			$L_{\zeta} \Rightarrow$	M ₁	L ₃			

**Radiative Partition Function (Γ) Governs the Relative Intensities
Nominal Values (Varies slowly with Atomic Number)**

K Shell	L Shell	M Shell
$K_{\alpha 1} = 100$	$L_{\alpha 1} = 100$	$M_{\alpha 1\gamma} = 100$
$K_{\alpha 2} = 50$	$L_{\alpha 2} = 50$	$M_{\beta} = 60$
$K_{\beta 1} = 15-30$	$L_{\beta 1} = 50$	
$K_{\beta 2} = 1-10$	$L_{\beta 2} = 20$	
$K_{\beta 3} = 6-15$	$L_{\beta 3} = 1-6$	
	$L_{\beta 4} = 3-5$	
	$L_{\gamma 1} = 1-10$	
	$L_{\gamma 3} = 0.5-2$	
	$L_{\eta} = 1$	
	$L_{\chi} = 1-3$	

Special Considerations for L Intensity Calculations:

$$\begin{aligned}
 I(L\alpha) &= \Gamma_{L3}(\alpha) \cdot \omega_{L3} \cdot Q_{L3} \\
 I(L\beta) &= \Gamma_{L3}(\beta 2) \cdot \omega_{L3} \cdot Q_{L3} + \Gamma_{L2}(\beta 1) \cdot \omega_{L2} \cdot Q_{L2} + \Gamma_{L1}(\beta 34) \cdot \omega_{L1} \cdot Q_{L1} \\
 I(L\gamma) &= \Gamma_{L3}(\gamma) \cdot \omega_{L3} \cdot Q_{L3} \\
 I(L\eta) &= \Gamma_{L2}(\eta) \cdot \omega_{L2} \cdot Q_{L2} \\
 I(L\chi) &= \Gamma_{L2}(\chi 1) \cdot \omega_{L2} \cdot Q_{L2} + \Gamma_{L1}(\chi 3) \cdot \omega_{L1} \cdot Q_{L1}
 \end{aligned}$$

The cross-sections used should also include

Coster-Kronig Transitions (F) and K Shell Vacancies (N)

$$Q_{L3}^{Total} = Q_{L3} + F_{L23} \cdot Q_{L2} + (F_{L13} + F_{L12} \cdot F_{L23}) \cdot Q_{L1} + N_{KL3} \cdot Q_K$$

$$Q_{L2}^{Total} = Q_{L2} + F_{L12} \cdot Q_{L1} + N_{KL2} \cdot Q_K$$

$$Q_{L1}^{Total} = Q_{L1} + F_{KL1} \cdot Q_K$$

Coster-Kronig Transitions:

This is an inverse transition where electrons travel up the potential well from a lower energy state to a higher one.

K Shell Vacancies:

Generation of a K-shell X-ray results from L shell electron dropping down to the K shell and an L-shell Vacancy is created. This will lead to a indirect emission of an L shell X-ray even though the L shell core hole was not initially created by direct electron excitation.

Quantitative Analysis using XEDS Standardless Method

Invoke the Intensity Ratio Method, that is consider the ratio of x-ray lines from two



$$\frac{I_A}{I_B} = \frac{\kappa_A \epsilon_A C_A}{\kappa_B \epsilon_B C_B} = k_{AB}^{-1} \frac{C_A}{C_B}$$

$$\kappa_A = \frac{\sigma_A \omega_A \Gamma_A}{W_A}$$

$$\frac{\kappa_A \epsilon_A}{\kappa_B \epsilon_B} = k_{AB}^{-1} \quad (\text{k-factor})$$

This simple equation states that the relative intensity ratio of any two characteristic x-ray lines is directly proportional to the relative composition ratio of their elemental components multiplied by some "constants" and is independent of thickness.

NOTE: The k_{AB} factor is not a universal constant!!
Only the ratio of κ_A/κ_B is a true physical constant and is independent of the AEM system. The ratio of ϵ_A/ϵ_B is not a constant since no two detectors are identical over their entire operational range. This can cause problems in some cases as we shall see.



The analysis to this point has only yielded the **relative compositions** of the specimen. We need one additional assumption to convert the relative intensity ratio's (I_i/I_j) into compositions namely:

$$\sum_{i=1}^N C_i = 1.0$$

One now has a set of N equations and N unknowns which be solved algebraically solved for the individual composition values.

Thus for a simple two element system we have:

$$\frac{I_A}{I_B} = k_{AB} \frac{C_A}{C_B}$$

and


$$C_A + C_B = 1.0$$

or

$$C_B \left(\frac{C_A}{C_B} + 1 \right) = 1$$

Solving for C_B and C_A

$$C_B = \left(\frac{1}{1 + \frac{C_A}{C_B}} \right) = \left(\frac{1}{1 + \frac{I_A}{I_B} \cdot k_{AB}} \right) \text{ and } C_A = C_B - 1$$




Vitek et al, AEM, 1984
Variation in Measured Composition on 308 SS for Different Labs

Example in which K-factor is stable

Cr, Fe, Ni

Note: Detector efficiency ~ 100% in this energy range

Variation in K-factor with AEM/Detector System Specimen: Uniform NiO film on Be Grid




Instrument	Experimental K-Factor	Apparent Variation in Composition
1	5.17	80.9/19.1
2	1.56	56.1/43.9
3	1.25	50.6/49.4
4	1.25	50.6/49.4
5	1.24	50.4/49.5
6	1.13	48.1/51.9
7	1.10	47.4/52.6
8	1.03	45.8/54.2
9	0.74	37.8/62.2

From: Comparison of UTW/WL X-ray Detectors on TEM/STEMs and STEMs

Thomas, Charlot, Franti, Garratt-Reed, Goodhew, Joy, Lee, Ng, Plieta, Zaluzec.
Analytical Electron Microscopy-1984

Variation in K-factor with AEM/Detector System

Low Energy End will not be the only problematic area



Determining the k_{AB}^{-1} Factor

Experimental Measurements

Prepare thin-film standards of known composition
then measure relative intensities and solve explicitly
for the k_{AB} factor needed. Prepare a working data base.

This is the "best" method, but

- specimen composition must be verified independently
- must have a standard for every element to be studied

Theoretical Calculations


Attempt first principles calculation knowing
some fundamental parameters of the AEM system



Start with a limited number of k_{AB} factor measurements,
then fit the AEM parameters to best match the data.
Extrapolate to systems where measurements and/or
standards do not exist.





Method 1. (Goldstein et al) Assume values for ρ , μ , and determine the best s to fit k_{AB} . This procedure essentially iterates the fit of s to the data.

Method 2. (Zaluzec) Assume values for ρ , μ , determine the best e to fit k_{AB} . This procedure essentially iterates the fit of e (detector window parameters) to the data.

K-Factor Calculation Experiment vs Theory







Sources of values for k_{AB} Calculations

W	- International Tables of Atomic Weights
$\Gamma(K)$	- Schreiber and Wims, X-ray Spectroscopy (1982) Vol 11, p. 42
$\Gamma(L)$	- Scofield, Atomic and Nuclear Data Tables (1974) Vol 14, #2, p. 121
$\omega(K)$	- Bambynek et al, Rev. Mod. Physics, Vol 44, p. 716 Freund, X-ray Spectrometry, (1975) Vol 4, p.90
$\omega(L)$	- Krause, J. Phys. Chem. Ref. Data (1974) Vol 8, p.307
$\sigma(E_0)$	- Inokuti, Rev. Mod. Physics, 43, No. 3, 297 (1971) - Goldstein et al, SEM 1, 315, (1977) - Chapman et al, X-ray Spectrometry, 12,153,(1983) - Rez, X-ray Spectrometry, 13, 55, (1984) - Egerton, Ultramicroscopy, 4, 169, (1969) - Zaluzec, AEM-1984, San Fran. Press. 279, (1984)
$\epsilon(E)$	- Use mass absorption coefficients from: -Thinh and Leroux; X-ray Spect. (1979), 8, p. 963 -Henke and Ebsiu, Adv. in X-ray Analysis,17, (1974) -Holton and Zaluzec, AEM-1984, San Fran Press,353,(1984)

Quantitative Analysis using XEDS Thin Film Standards Method

Invoke the Intensity Ratio Method, but now consider the ratio of the same x-ray line from two different specimens, where one is from a **standard** of known composition while the other is **unknown**:

$$\frac{I_u}{I_s} = \frac{\eta_u \rho_u t_u}{\eta_s \rho_s t_s} * \frac{C_u}{C_s}$$

$$C_u = \frac{\eta_s \rho_s t_s}{\eta_u \rho_u t_u} * \frac{I_u}{I_s} * C_s$$

This simple equation states that the relative intensity ratio of same characteristic x-ray line is directly proportional to the relative composition ratio of the two specimens multiplied by a some new parameters.

η = incident beam current
 ρ = local specimen density
 t = local specimen thickness

Quantitative Analysis using XEDS Specimen Thickness Effects

For finite thickness specimens, what is a thin film?

Previous Assumptions:

- No Energy loss,
- No X-ray absorption,
- No X-ray fluorescence

NOTE: Electron Transparency is insufficient!

Effects of energy loss on Characteristic X-ray Production:

$$I_A = \frac{N_O \rho C_A}{W_A} \omega_A \Gamma_A \eta_p \Omega \epsilon_A \int_0^{t_o} \sigma_A(E) dt$$

$$I_A = \frac{N_O \rho C_A}{W_A} \omega_A \Gamma_A \eta_p \Omega \epsilon_A \int_0^{t_o} \sigma_A(E) dt$$

Calculated Change in σ with Energy Loss

$$\frac{\Delta\sigma}{\sigma} = \frac{\sigma(100) - \sigma(100 - \Delta E)}{\sigma(100)}$$

Energy Loss(eV)	Al	Ni	Ag
500	0.004	0.003	0.002
1000	0.009	0.006	0.004
2000	0.018	0.013	0.009
5000	0.047	0.030	0.022

Quantitative Analysis using XEDS Specimen Thickness Effects

Consider the mean energy loss of a electron beam () going through a film of thickness = "t". Let s = electron pathlength (~thickness) then (Inokuti -1971):

$$\frac{dE}{ds} = \frac{ne^4 N_0 \rho Z}{WE} \left\{ \ln \left(\frac{2ET_E}{J^2(1-\beta^2)} \right) - 2 \ln(2) \sqrt{(1-\beta^2)} + 1 - \beta^2 + \frac{(1-\sqrt{1-\beta^2})^2}{8} \right\}$$


Calculated Energy Losses (eV) as a function of pathlength

Voltage	Pathlength	Al	Ni
100 kV	100 nm	87	255
	400 nm	349	1020
300 kV	100 nm	51	150
	400 nm	203	60

For conventional specimen thicknesses used in AEM, the effect of energy loss on characteristic x-ray production is negligible.

Quantitative Analysis using XEDS : Absorption Correction

$$I = \int_0^{t_0} I_0(t) \exp(-\mu^* d) dt = \int_0^{t_0} I_0(t) \exp(-\chi \rho t) dt$$




$$d = \text{Absorption Pathlength from the specimen to detector} = t \frac{\sin(\beta)}{\cos(\beta - \theta_E)}$$

χ = Geometrical factor multiplied by average mass absorption coefficient for the measured x-ray line in the compound

$$= \left(\frac{\mu}{\rho} \right)_{\text{Spec.}} \frac{\sin(\beta)}{\cos(\beta - \theta_E)}$$

$$\left(\frac{\mu}{\rho} \right)_{\text{Spec.}} = \text{Weighted average mass absorption coefficient}$$

$$= \sum_{i=1}^N \left(\frac{\mu}{\rho} \right)_i \cdot C_i \text{ (Note: composition dependent!!)}$$



$$\beta = \text{Arccos} \left(\frac{\sin\phi_y \cdot \cos\theta_A \cdot \cos\theta_E + \sin\phi_x \cdot \cos\phi_y \cdot \cos\theta_E \cdot \sin\theta_A}{\sqrt{a^2 + b^2 + c^2}} \right)$$

$$a = \cos\phi_x \cdot \cos\phi_y \cdot \cos\theta_A \cdot \cos\theta_E$$

$$b = \cos\phi_x \cdot \cos\phi_y \cdot \cos\theta_E \cdot \sin\theta_A$$

$$c = \sin\phi_y \cdot \cos\theta_A \cdot \cos\theta_E + \sin\phi_x \cdot \cos\phi_y \cdot \cos\theta_E \cdot \sin\theta_A$$

Now rederive the standardless equations to include absorption.

$$\frac{I_A}{I_B} = \frac{\epsilon_A}{\epsilon_B} * \frac{\kappa_A}{\kappa_B} * \frac{\delta_A}{\delta_B} * \frac{C_A}{C_B}$$

with

$$\frac{\delta_A}{\delta_B} = \frac{\left(\frac{\mu}{\rho}\right)_{AB}^{B \text{ in}} (1 - \exp(-\chi \rho t^*)_{AB}^{A \text{ in}})}{\left(\frac{\mu}{\rho}\right)_{AB}^{A \text{ in}} (1 - \exp(-\chi \rho t^*)_{AB}^{B \text{ in}})}$$

$$t^* = t_0 \frac{\sin(\beta)}{\cos(\beta - \theta_E)}$$

$$\chi = \left(\frac{\mu}{\rho}\right)_{\text{X-ray in Compound}} = \sum_{i=1}^N \left(\frac{\mu}{\rho}\right)_i * C_i$$


β = Electron Incidence Angle

= Function of Stage Tilts: ϕ_x, ϕ_y , & Detector Azimuth θ_A

θ_E = Detector Elevation Angle


Define the Thin Film approximation: $\chi \rho t^* < 0.1$

Thin Film approximation: $\chi \rho t^* < 0.1$



$$t^* \leq \frac{0.1}{\chi \rho}$$

Effect of Tilting Specimen on Absorption



For a plane parallel slab specimen, tilting has the effect of increasing the Specimen Thickness.
Different Detector /Specimen Geometries will enhance/reduce the Absorption Effects




FIG. 1.--Theoretically calculated intensity ratio of I_0^2/I_{00} as function of specimen tilt ϕ_A at detector/specimen geometry of $\theta_E = 0^\circ$, $\theta_A = 45^\circ$ for specimen thickness values ranging from 10 to 400 nm as indicated.




FIG. 2.--Theoretically calculated intensity ratio of I_0^2/I_{00} as function of specimen tilt ϕ_A at detector/specimen geometry of $\theta_E = 20^\circ$, $\theta_A = 30^\circ$ for specimen thickness values ranging from 10 to 400 nm as indicated.






FIG. 3.--Theoretically calculated intensity ratio of I_0^2/I_{00} as function of specimen tilt ϕ_A at detector/specimen geometry of $\theta_E = 60^\circ$, $\theta_A = 0^\circ$ for specimen thickness values ranging from 10 to 400 nm as indicated.

In the "real" world few specimens have the shape used to derive this correction. Next consider two representative geometries:

Symmetric Wedge Geometry:



Replace all β 's by $\beta \pm \alpha/2$, where α is the wedge angle of the specimen.
 $+\alpha/2$ variant applies when the detector is positioned such that the pathlength *increases* relative to the parallel slab model.
 $-\alpha/2$ variant applies when the detector is positioned such that the pathlength *decreases* relative to the parallel slab model.



Glitz et al (MAS-1981)

Attempt a Wedge Model Correction using previous formulae.

Path #1		Path #2	
$\frac{I_{Ni}}{I_{Al}}$	= 2.64	$\frac{I_{Ni}}{I_{Al}}$	= 5.18
Thin Film Model	$\frac{Ni=60.9}{Al=39.1}$	Thin Film Model	Error!
Parallel Slab	$\frac{Ni=45.3}{Al=54.7}$	Parallel Slab	$\frac{Ni=60.5}{Al=39.5}$
Wedge Model	$\frac{Ni=50.6}{Al=49.4}$	Wedge Model	$\frac{Ni=42.6}{Al=57.4}$
$\chi_{pt} Ni$	= 0.016	$\chi_{pt} Ni$	= 0.081
$\chi_{pt} Al$	= 0.925	$\chi_{pt} Al$	= 4.24

Absorption Correction has limited applications keep $\chi_{pt} < 1$

Specimen Homogeneity


In this and all other derivations we have assumed that over the excited volume, as well as along the exiting pathlength, the specimen is homogeneous in composition. If this assumption is invalid, one must reformulate the absorption correction and take into account changes in μ/ρ , Z , and t along the exiting pathlength.

Effects of Beam Broadening


Parallel Slab Model: No Change in absorption pathlength
 Wedge Model: There is a correction the magnitude of which varies with the wedge angle.

Effects of Irregular Surface

This cannot be analytically modeled but must be understood!



X-Ray Fluorescence Correction



$$I_A^{\text{Measured}} = I_A^{\text{Electron}} + I_A^{\text{XRF by B}}$$

$$I_A^{\text{Measured}} = I_A^{\text{Electron}} \left(1 + \frac{I_A^{\text{XRF by B}}}{I_A^{\text{Electron}}} \right)$$

Calculation of the Intensity of I_A x-rays due to fluorescence by I_B

$$dI_A^{XRF \text{ by } B} = I_B \cdot \left[\exp \left(-\left(\frac{\mu^B}{\rho}\right)_{AB} \rho z \operatorname{cosec} \phi \right) \right] \cdot \left[\rho C_A \left(\frac{\mu^B}{\rho}\right)_A \operatorname{cosec}(\phi) \right] \cdot \left[\left(\frac{\Delta\mu}{\mu}\right)_A \right] \cdot \left[\omega_A \cdot \Gamma_A \right]$$

$$dI_A^{XRF \text{ by } B} = I_B \cdot P_A^{XRF \text{ by } B} dz d\phi$$

$$\left[\exp \left(-\left(\frac{\mu^B}{\rho}\right)_{AB} \rho z \operatorname{cosec} \phi \right) \right]$$

Fraction of B x-rays that are absorbed in the material AB on leaving at angle ϕ from depth z

$$\left[\rho C_A \left(\frac{\mu^B}{\rho}\right)_A \operatorname{cosec}(\phi) \right]$$

Probability that B x-rays are absorbed by atoms of type A

$$\left[\left(\frac{\Delta\mu}{\mu}\right)_A \right]$$

Probability that upon absorption the K-Shell is ionized for A

$$\left[\omega_A \cdot \Gamma_A \right]$$

Probability that a $K\alpha$ x-ray is emitted

Now integrate this over the specimen geometry $\phi=4\pi$, T =thickness

$$I_A^{XRF} = \int_0^T \int_0^{4\pi} dI_A^{XRF} = \int_0^T \int_0^{4\pi} I_B \cdot P_A^{XRF} dz d\phi$$

Substituting for I_B and dividing by I_A and generalizing to element i

$$\frac{I_A^{XRF \text{ by } i}}{I_A} = C_i \omega_i \Gamma_i \frac{A_A}{A_i} \frac{Q_i}{Q_A} \left(\frac{\Delta\mu}{\mu}\right)_A \left(\frac{\mu^i}{\rho}\right)_{Ai} \int_0^T \int_0^{4\pi} (\text{terms in Exp. \& cosec } \phi) dz d\phi$$

Solution depends on the geometric assumptions
(see for example Twigg & Fraser MAS, 1982)

$$\frac{I_A^{XRF \text{ by } i}}{I_A} = C_i \omega_i \Gamma_i \frac{A_A}{A_i} \frac{Q_i}{Q_A} \left(\frac{\Delta\mu}{\mu}\right)_A \left(\frac{\mu^i}{\rho}\right)_A \frac{\rho T}{2} \left[1.12 + \left(\frac{\mu^i}{\rho}\right)_{Ai} \cdot \frac{\rho T}{4} - \ln \left(\left(\frac{\mu^i}{\rho}\right)_{Ai} \cdot \rho T \right) \right]$$

for a plane parallel slab specimen

Now

$$\text{Define } \gamma_A = \left(1 + \sum_{i=1}^N \frac{I_A^{XRF \text{ by } i}}{I_A^{\text{Electron}}} \right)$$

Next, rederive the standardless equations to include x-ray fluorescence and you can show:

$$\frac{I_A}{I_B} = \frac{\epsilon_A}{\epsilon_B} \times \frac{\kappa_A}{\kappa_B} \times \frac{\delta_A}{\delta_B} \times \frac{\gamma_A}{\gamma_B} \times \frac{C_A}{C_B}$$

as in the case of x-ray absorption this requires iterative solution because the ratio of γ 's are composition dependent.

When is the XRF Correction important?

- When fluorescing line is near the absorption edge of the lower energy line. Typically within a few atomic numbers (i.e. Z+2 to Z+6)
- When specimen is thick or path length is long

Define a thin film approximation for XRF as:

$$\frac{I_A^{XRF \text{ by } i}}{I_A^{Electron}} < 0.05$$

$$C_i \propto \Gamma_i \frac{A_A}{A_i} \frac{Q_i}{Q_A} \left(\frac{\Delta\mu}{\mu} \right)_A \left(\frac{\mu^i}{\rho} \right)_A \frac{\rho T}{2} \left[1.12 + \left(\frac{\mu^i}{\rho} \right)_A \cdot \frac{\rho T}{4} - \ln \left(\left(\frac{\mu^i}{\rho} \right)_A \cdot \rho T \right) \right] < 0.05$$

Examples Absorption & Fluorescence Corrections

Assume 100 kV, $\theta_A = 90^\circ$, $\theta_E = 20^\circ$, $\beta = 35^\circ$, $\alpha = 0^\circ$ i.e. tilted parallel slab

Ni 90%, Fe 10%

Calculation of Thin Film approximation $\chi_{pt} \Rightarrow 1.25\mu\text{m}$ for Fe in the Alloy

\therefore for TEM specimens we can almost always ignore absorption effects

What about XRF?

$$\text{Let } T = 5000\text{\AA} \quad \frac{I_{Fe}^{XRF \text{ by Ni}}}{I_{Fe}^{Electron}} = 0.103 \text{ i.e. } \sim 10\%$$

$$\text{Let } T = 1000\text{\AA} \quad \frac{I_{Fe}^{XRF \text{ by Ni}}}{I_{Fe}^{Electron}} = 0.028 \text{ i.e. } 2.8\%$$

\therefore for TEM specimens we may be affected by XRF

Ni 90%, Al 10%

Calculation of Thin Film approximation $\chi_{pt} \Rightarrow 232 \text{ \AA}$ for Al in Alloy!!!

\therefore for TEM specimens we cannot ignore absorption effects

What about XRF?

Let $T = 5000 \text{ \AA}$

$$\frac{I_{Al}^{XRF \text{ by Ni}}}{I_{Al}^{Electron}} = 0.000186 \text{ i.e. } 0.01\%$$

\therefore for typical TEM specimens we can ignore XRF effects

Quantitative Analysis using XEDS - Bulk

For a Bulk Specimen:

$$I_A^{K\alpha} = \left[\Gamma_A \omega_A \int_{\epsilon}^{\epsilon} \frac{\sigma_A(E, Z)}{S_A} dE \right] R_A \{f(\chi)\} \{F_A\} C_A \left[\frac{N_0 \rho}{W_A} \right] \eta_0 \epsilon_A \Omega$$

- I_A = Measured x-ray intensity per unit area
- σ = K^{th} -shell ionization cross-section
- ω = K^{th} -shell fluorescence yield
- Γ = K^{th} -shell radiative partition function
- R = Backscatter Correction Term
- S = Electron Stopping Power
- $f(\chi)$ = Absorption Correction Term
- F = Fluorescence Correction Term
- W = Atomic Weight
- N_0 = Avagadro's number
- ρ = Density
- C = Composition (At %)
- η_0 = Incident electron flux
- t = Specimen thickness
- ϵ = Detector efficiency
- Ω = Detector solid angle

Quantitative Analysis using XEDS Standards Method

Invoke the Intensity Ratio Method, but now consider the ratio of the same x-ray line from two different specimens, where one is from a **standard** of known composition while the other is **unknown**:

$$\frac{I_U^{K\alpha}}{I_S^{K\alpha}} = \frac{\left\{ \Gamma_A \omega_A \int_{E_c}^E \frac{\sigma_A(E,Z)}{S_U} dE \right\} \{R_U\} \{f_U(\lambda)\} \{F_U\} C_U \left\{ \frac{N_0 \rho}{W_A} \right\} \{\tau_0 \epsilon_A \Omega\}}{\left\{ \Gamma_A \omega_A \int_{E_c}^E \frac{\sigma_A(E,Z)}{S_S} dE \right\} \{R_S\} \{f_S(\lambda)\} \{F_S\} C_S \left\{ \frac{N_0 \rho}{W_A} \right\} \{\tau_0 \epsilon_A \Omega\}}$$

$$\frac{I_U^{K\alpha}}{I_S^{K\alpha}} \equiv K_{\text{Ratio}} = K_Z \cdot K_A \cdot K_F \cdot \frac{C_U}{C_S}$$

or rearranging

$$C_U = K_{\text{Ratio}} \cdot (\text{ZAF Corrections}) \cdot C_S$$

These equations state that the relative intensity ratio of same characteristic x-ray line is directly proportional to the relative composition ratio of the two specimens multiplied by a some correction terms

Advantages

Standard may be a pure element.
The closer the standard is to the unknown material the smaller the correction and higher the accuracy.

Disadvantages


Effects of surface films can be critical
Must have a standard for each element to be analyzed


Additional Topics

Heterogeneous Specimens
Composition Profiles
Electron Channeling
Radiation Damage

All quantitative analysis equations were derived assuming that the specimen is homogeneous over the excited volume

Application of these equations to heterogeneous specimens effectively averages the composition over the excited volume.





$$C^*(x,y) = C(x,y,z) * d(x,y,z)$$

$C^*(x,y)$ = Apparent profile measured

$C(x,y,z)$ = Actual composition profile

$d(x,y,z)$ = Incident beam profile

$*$ = Convolution operator


F, F^{-1} = Fourier and Inverse Fourier Transforms

In the 2 dimensional limit one can deconvolute the measured profile using:


$$C(x,y) = F^{-1} \left[\frac{F[C^*(x,y)]}{F[d(x,y)]} \right]$$

Realistically, it is better to decrease the probe diameter and specimen thickness


Electron Channeling Effects





Orientation Dependence of Intensity Ratio in Ordered System



High Angular Resolution Electron Channeling X-ray Spectroscopy




ALCHEMI Atom Location by CHanneling EMission Tafto & Spence - Science 1982




$(\text{Mg}_{0.90}\text{Fe}_{0.10}\text{Ni}_{0.004}\text{Mn}_{0.002})_2\text{SiO}_4$
Note: Mn prefers Si Sites and Ni prefers Mg Sites

Studies of Spinel



with T. Soeda, S. Matsumura, L. Rehn
Kyushu Univ. / ANL




MgAl_2O_4 Cross-Section


HAREC XS -



Mg
Al
O




{110}



with T. Soeda, S. Matsumura, L. Rehn
Kyushu Univ. / ANL

Coherent Bremsstrahlung
Reese, Spence, and Yamamoto- Phil Mag. 1984 Vol 49, 5, 697



Radiation Damage vs Microanalysis

Comparison of Characteristic Signal Generation and Displacement Rates

$$\frac{\text{Ionizations}}{\text{sec}} \quad N_i = Q * J * C \frac{N_{op}}{A}$$

$$\frac{\text{Displacements}}{\text{sec}} \quad N_d = \sigma_d * J * \phi(f) * C \frac{N_{op}}{A}$$

$$\frac{\text{Detected X-rays}}{\text{sec}} \quad N_i = Q * [\omega_X \Gamma_X \epsilon_X \Omega_X] * J * C \frac{N_{op}}{A}$$

$$\frac{\text{Detected EEL}}{\text{sec}} \quad N_i = Q * [\omega_E \Gamma_E \epsilon_E \Omega_E] * J * C \frac{N_{op}}{A}$$

$$\omega_E \sim \Gamma_E \sim \epsilon_E \sim 1$$

Radiation Damage vs Microanalysis

Energy Transferred by an Incident Electron to an Atomic Nucleus

$$T_T = \frac{2 \cdot T_0 \cdot (T_0 + 2 \cdot m_0 c^2) \cdot \sin^2(\frac{\theta}{2})}{M c^2}$$

$T_0 = eV_0$, $M =$ nuclear mass, $m_0 =$ electron mass, $\theta =$ scattering angle

Comparison of Maximum Transferable Kinetic Energy for Selected Elements with Displacement and Sputtering Energies (all values in eV)

Element	100 kV	200 kV	300 kV	400 kV	T_d^1	T_s^2
Al	8.93	19.5	31.6	45.3	16	3.5-7.0
Ti	5.00	11.0	17.8	25.5	15	4.9-9.8
V	4.73	10.3	16.72	24.0	29	5.3-10.6
Cr	4.63	10.1	16.38	23.5	21	4.1-8.2
Fe	4.31	9.40	15.25	21.8	16	4.3-8.6
Co	4.08	8.91	14.45	20.7	23	4.4-8.8
Ni	4.10	8.94	14.5	20.8	21	4.5-9.0
Cu	3.79	8.26	13.4	19.2	18	3.5-7.0
Zn	3.69	8.03	13.03	18.7	16	1.4-2.8
Nb	2.59	5.65	9.17	13.2	24	7.5-15.0
Mo	2.51	5.47	8.88	12.7	27	6.8-13.6
Ag	2.23	4.87	7.90	11.3	28	3.0-6.0
Cd	2.14	4.67	7.58	10.9	20	1.2-2.4
Ta	1.33	2.90	4.71	6.75	33	8.1-16.2
Pt	1.23	2.69	4.37	6.26	33	5.9-11.8
Au	1.22	2.67	4.32	6.2	36	3.8-7.6

T_d = Displacement Threshold (Bulk) Energy

T_s = Sputtering Threshold (Surface) Energy

Calculated Sputtering Cross-Sections

